

# **PREPARATION OF DEFECT-FREE POLYACRYLAMIDE ELECTROPHORESIS GELS IN PLASTIC CASSETTES**

## **ABSTRACT OF THE DISCLOSURE**

Irregularities in the pore structures of polyacrylamide gels that are formed in cassettes against plastic walls are reduced or eliminated by the inclusion of an oxygen scavenger in the gel-forming solution. Avoidance of the irregularities results in electropherograms with fewer distortions in the solute bands.

SF 1301098 v2

10087140-022702